

SafeSense[®]



Bed Exit assistance system

SafeSense® for greater freedom

This innovative early-warning system to notify staff when the patient leaves the bed is a reliable aid in everyday nursing care. The smart bed sensor technology helps to prevent falls and ensure greater mobility for the patient.




24/7 patient control

SAFE PATIENT MOBILITY

The nursing staff can immediately respond and assist the patient.

After the set interval has elapsed, the SafeSense® system will send a signal to call staff/nurses.

- Prompt assistance if the patient falls.
- A high degree of safety during night and weekend shifts, when there are fewer members of staff.
- The system can also trigger the alarm if the patient tries to climb over the side rails.




Falls are the most common cause of injury or death for people over the age of 65*

and safety

FOR BETTER ORIENTATION

When the patient gets out of the bed, the built-in night light automatically turns on, which greatly facilitates orientation and makes the patient more independent. After the patient returns to the bed, the light switches off automatically.



FLEXIBLE TIMING


The timing of the system for calling staff can be set as an instant alarm (0 seconds) or at intervals of up to 30 minutes.

- Flexible use after the patient's physical and mental condition have been assessed.
- Individual settings to suit the patient's degree of mobility.



SYSTEM INTEGRATION

The system is quick and easy to install using P&P connectors on a compatible bed, as needed.



SafeSense®

Patient

- Promotes independent mobility
- Minimises the risk of falls
- Automatic night light



- Staff are quickly informed if the patient leaves the bed
- Option to use individual settings
- Easy to control



Staff

Management

- A high degree of safety
- Compatible with system for calling nurses
- Minimal maintenance costs



Members of LINET Group

LINET spol. s r.o.

Želevčice 5 | 274 01 Slaný | Czech Republic
tel.: +420 312 576 400 | fax: +420 312 522 668 | e-mail: info@linet.com | www.linet.com

