


■ ■ wissner-
■ ■ bosserhoff



Accessories

Nursing care





Accessories

With love for detail

When it comes to care for the elderly, a good human interaction between nursing staff and residents is very important. wissner-bosserhoff understood that an open-minded and friendly personal interaction can be supported by helpful and innovative products. When less time is needed for necessary jobs and tasks, more time is left to care for residents. We have developed useful accessories for our products, which contribute to better quality and efficiency in the field of nursing care.

Contents

Lights	4 - 5
Bed exit system	6 - 7
Mobilization	8 - 11
Storage	12 - 13
Protection	14 - 17
Care	18 - 19
Design Softcover	20 - 23



Lights


The undulating shape of the lampshade, with its partially translucent character, is prominent as a special design element of soluna lights. Equipped with a swivel arm, the light allows the light focus to be precisely positioned close to the object to be illuminated, regardless of whether the light is attached to the bed or the bedside cabinet. Either halogen or LED illuminants can be chosen as the light source. Advantage of LEDs: a long life, low energy consumption, no heat generated and absolutely safe operation.



50-0677
LED reading light lymera
(with and without power supply unit)
Transparent lampshade



50-0789
Design reading light soluna
HALOGEN, incl. 12 V transformer
50-0550
Design reading light soluna
LED, without power supply unit



01-010812, 01-010813
Switchable underbed light
(Night light function)

Advantages of LED

- Long service life
- Low energy consumption
- No heat build-up
- Absolutely safe operation



Safety by day and night

SafeSense® – the intelligent bed exit sensor system



Video SafeSense®

wissner-bosserhoff, your reliable partner for the reduction of freedom-depriving measures (FDM): With SafeSense®, wissner-bosserhoff presents an intelligent assistance sensor system, which reliably helps reducing freedom-depriving measures day and night. It is equipped with an automatic night light, which is activated when the resident leaves the bed and is deactivated when the resident is back in bed again. SafeSense® is compatible with standard nurse call systems and can be retrofitted to wissner-bosserhoff nursing home beds manufactured in 2002 or later.



Automatic night light

- Night light is automatically switched on when the resident leaves the bed
- Night light is switched off again as soon as the resident is back in bed

Individual alarm system

- Freely adjustable immediate alarm from 0 seconds to 30 minutes
- Adaptation to the resident's mobility

Integrated or available as retrofit kit*

- Accessory option for current WIBO nursing home beds
- Retrofit kit for electrical WIBO nursing home beds manufactured in 2002 or later
- Available with cable or wireless




More safety –
less freedom-depriving
measures (FDM)


- Faster response time in the event of falls
- Greater safety during night and weekend shifts
- Fewer inspections
- Ideal in combination with sentida low nursing beds and SafeFree® side guards

SafeSense® bed exit system

SafeSense® Sensor mat 2.1




02-000876
SafeSense® sensor mat 2.1
80 x 20 cm



04-1851
Digital bed exit box 2
cable or wireless version

**Also available
as a set:**
Sensor Mat 2.1 &
Bed Exit Box 2
10-0871
Cable version
10-0872
Radio version

SafeSense® Connection kit




04-1508
Adapter cable with on/off switch and suitable
plug for connection to the nurse call system
(customer-specific) for the cable and
wireless version




04-1608
Adapter cable without on/off switch
with suitable plug to connect to the
nurse call system (customer-specific)
for the cable and wireless version

SafeSense® Retaining plate



01-011193
Retaining plate for DU09 motor



01-011194
Retaining plate for Q7 motor

Mobilizing


Independent mobilization reduces the risk of falling. Our accessories help the resident to maintain their independence. We would all like to remain independently mobile for as long as possible. This is why it is so important to preserve and promote the individual mobility of care-dependent residents in line with the standard, developed by experts, known as „Preservation and Promotion of Mobility in the Field of Care“. The benefits for resident range from social participation to health promotion.



SafeLift bed operation




10-0811
SafeLift Controller




10-0808
SafeLift Box

As a set:
10-0799 SafeLift set

SafeLift Retaining plate



01-011193
for DU09 motor



01-011194
for Q7 motor



10-0799
SafeLift Mobilisations-Assistent,
Set including box

- Maintain & improve mobility
- Simple bed operation
- No cable, no battery
- An automatic night light
- Easy installation, retrofitting

Mehr erfahren:



The MobiStick 2 enables clever bed operation for residents. Its ergonomic shape makes it easy and intuitive to use anywhere. It encourages individual mobilization. Residents can effectively support and move themselves.



MobiStick 2

NEW!



Without tools

MobiStick 2 can be attached without tools on both sides of the bed, is suitable for every resident size and can be retrofitted for wissner-bosserhoff beds from 2004 onwards.



Due to its ergonomic shape MobiStick2 is easy and intuitive to use and universally applicable. The hand knob is rotatable and has a simple locking.



When not in use, it can simply be folded away and is therefore not in the way, but quickly ready for use when needed.

Learn more:



Mobilizing



10-0586

Serving tray with bookholding rails



10-0715

Standing-up aid for 09 and 01 side guards, as well as for ¾ side guards



10-0728

Accessory adapter for standing-up aid (10-0715)





04-0674 / 10-0502

Lifting pole with oval trapeze handle and belt retractor



Quick and convenient access
for personnel and residents



10-0422 (A)

Accessory holder to attach to continuous
(09) side guards

10-0421 (B)

Handset holder can be plugged into the
side guards accessory holder (10-0422)



10-0419 (A)

Adapter to attach to the divided (10)
side guards

10-0421 (B)

Handset holder can be plugged into the
side guards accessory holder (10-0419)

Storage

An extensive range of accessories to supplement our products and to support everyday nursing care completes the overall portfolio. This makes nursing care considerably easier for the staff and valuable time is saved, which can then be used to increase the quality of care.



10-0617
Integrated linen holder,
compatible with all beds




Integrated linen holder, pulled out



10-0719
Linen holder



10-0563/10-0533/10-0521
Segufix holder for thigh rest
Segufix holder for the backrest
Segufix holder for the seat section



Protector to close the foot-end gap of 220 cm beds

10-0580 for extension via cable pull

10-0570 for patient surface pull-out screwed tight

10-0628 for patient surface width of 100 cm and patient surface pull-out screwed tight



10-0595

Patient surface extension for beds with integrated extension, wooden decor



**Bed extension
for beds without integrated patient surface
extension**

10-0497

Extension element to extend the lower leg rest by 20 cm

10-0573

Extension element to extend the calf rest by 20 cm with 100 cm wide patient surface

Protection

Depending on how the bed is used, it is possible to prevent the resident from injuries and the bed from damage caused by the surroundings. Side guard cushioning compliments the bed and closes the gaps.



50-0206
Fall mat optimal for sentida beds



10-0597
Side guard extension
for continuous (09) side guard

10-0605
Side guard extension
for divided (10) side guards (no image)



50-0544
Evacuation sheet with belt and fixation straps



50-0155

Side section padding with push buttons
for divided (10) side guards



50-0174

Side section padding for divided (10)
side guards with safety strap, closed type



50-0207

Side guard net to close the gaps
of divided (10) side guards



50-0338

Side section padding for continuous (09)
two-part side guards for bed length of 200 cm

50-0339

Side section padding for continuous (09)
two-part side guards for bed length of 220 cm
(no image)

Protection



10-0682

3D wall stopper castor, reinforced type

10-0683

3D wall stopper castor for design A (no image)



10-0485

Reinforced wall stopper castors for divided (10) side guards as well as continuous 2-part and 3-part 09 side guards



10-0626

Vertical wall stopper castors



Vertical wall stopper castors, extended type

10-0507 extended by 22 cm


19-0494 extended by 12 cm

19-0480 extended by 9 cm



02-000617

Side panel buffer for divided (10) side guards



10-0530

Head-end wall spacer bracket carisma
Distance to wall approx. 17.4 cm



10-0548

Head-end wall spacer bracket movita
Distance to wall approx. 8 cm



10-0587

Head-end wall spacer bracket movita casa
Distance to wall approx. 11.5 cm



Head-end wall spacer, telescopic,
only for sentida 5-le and 6-le

10-0632 short, 1 x 20 mm

10-0622 medium, 4 x 20 mm

10-0612 long, 8 x 20 mm




Side-mounted wall spacer, telescopic,
only for sentida 5-le and 6-le

10-0633 short, 1 x 20 mm

10-0623 medium, 4 x 20 mm

10-0613 long, 8 x 20 mm



Head-end wall spacer, telescopic,
only for sentida 5 and 6

10-0801 short, 1 x 20 mm

10-0802 medium, 4 x 20 mm

10-0803 long, 8 x 20 mm

Head-end wall spacer,
only for sentida 5 and 6 (no image)

10-0707 short

10-0709 long



Side-mounted wall spacer, telescopic,
only for sentida 5 and 6

10-0804 short, 1 x 20 mm

10-0805 medium, 4 x 20 mm

10-0806 long, 8 x 20 mm

Side-mounted wall spacer,
only for sentida 5 and 6 (no image)

10-0708 short

10-0710 long

Care

The accessories for the accessory holder provide useful support for patient and staff.
The holder can be attached to the left and right side of the bed.



50-0693
Bent infusion stand with 4 hooks,
height-adjustable



50-0003
Drip clasp to attach to the lifting pole



50-0694
Straight infusion stand with 4 hooks
and drip vessel **01-006094**



50-0108
Name tag holder



10-0424
Accessories holder for divided (10) side guards, to attach to the side panel



50-0345
Shortened urine bottle basket (similar to picture)



50-0164
Urine bag holder (depicted here with urine bag and accessories holder)



10-0511
Bed transport trolley for assembled nursing home beds



Design | Softcover

With our new cover solution made of high-quality materials from renowned partners we offer an excellent surface material with outstanding properties.





sentida 7-i | Design Df | Decor Cottage Pine + fabric cover imitation leather silver grey

New
Combination-
possibilities
We advise you're
welcome.

Softcover | Synthetic leather



The special feel, special properties and expressive embossing give high-quality upholstery materials their individual character. In addition to unusual designs, surfaces are also distinguished by their ease of care. Imitation leather made of PVC plastic scores above all with durability and thus stands for genuine quality. Further information is available on request.



Silver gray (F6461730)



Pearl (F6461662)




Birch (F6461661)



Pebble (F6461653)



Cashmere (F6461737)



Smoke (F6461660)



Cherry (F6461665)



Corn (F6461706)



Lime (F6461657)



Royal (F6461710)

Properties of all materials shown here:

- Resistant to moisture and dirt
- Breathable
- Water-repellent
- Blood and urine resistant
- Disinfectant resistant
- Extremely hard-wearing
- Environmentally friendly and non-toxic
- Durable
- Flame retardant

(in accordance with the ignitability criteria in accordance with DIN EN 1021 Part 1 + 2 and BS 5852 Crib 5*)

Softcover | Fabrics



These upholstery fabrics are unrivalled in their clean and neat appearance. The patented innovative system leaves bacteria, dirt and moisture no chance. Beautiful designs on first-class fabric qualities offer extensive combination possibilities between upholstery and decoration fabrics. Brand yarns made of Polyester FR or Trevira CS are used here. Further information is available on request.



Green, floral (Aurora 6520)



Blue, striped (Agon 5520)




Gray, striped (Agon 8520)



Gray, polka dot (Flint 8520)



Meadow (Mailo 6520)



Sea (Mailo 5520)



Sand (Mailo 1521)



Stone (Mailo 8520)



Rose (Mailo 4520)

50-0818

Softcover for the headboard Df for sentida healthcare beds

50-0819

Softcover for the foot section Df for beds with standard extension**

50-0820

Softcover for the foot section Df for beds with express extension**



Members of LINET Group

wissner-bosserhoff GmbH

Hauptstr. 4-6 | 58739 Wickede (Ruhr) | Germany

Phone: +49 2377.784-0 | Fax: -163 | info@wi-bo.de | wi-bo.com

